

IMPORTANT: PLEASE READ CAREFULLY AND KEEP THIS INFORMATION FOR FUTURE USE.

BAUSCH + LOMB

Boston[®]

One Step Liquid Enzymatic Cleaner

Specifically Formulated For Fluoro Silicone Acrylate and Silicone Acrylate Rigid Gas Permeable Contact Lenses

Boston[®] One Step Liquid Enzymatic Cleaner is a simple and convenient way to effectively remove protein deposits from rigid gas permeable contact lenses on a weekly basis. This clear, odorless enzymatic cleaner works together with the disinfecting step, right in your lens case, replacing enzymatic tablets and eliminating the need for separate vials. Use with either Boston[®] Conditioning Solution or Boston ADVANCE[®] Conditioning Solution.

DIRECTIONS (USE WEEKLY):

1. After daily cleaning (with Boston[®] Cleaner or Boston ADVANCE[®] Cleaner) and rinsing, place lenses in lens case and fill to top with any fresh Boston[®] Conditioning Solution or Boston ADVANCE[®] Comfort Formula Conditioning Solution. **It is important that the lens case is completely filled with fresh solution to ensure proper dilution of the Liquid Enzymatic Cleaner to reduce any chance of irritation.**
2. Detach one container from the strip by pulling tab toward you from top to bottom. To open, twist the tab off the single-use container. Squeeze container to release 2 drops of Boston[®] One Step Liquid Enzymatic Cleaner into each lens well already filled with Conditioning Solution. **Soak lenses for at least 4 hours (or overnight) before wearing and to**

ensure proper cleaning and disinfecting. Always use fresh solution for **soaking and storing lenses.** Discard container immediately after use.

3. After removing lenses from the lens case, rub and rinse lenses for approximately 5 seconds with any Boston[®] Conditioning Solution to ensure that the enzymatic cleaner is thoroughly rinsed off the lenses. **DO NOT USE TAP WATER FOR THIS FINAL RINSE.**
4. Apply fresh Conditioning Solution to wet the lenses for additional cushioning, if desired, and insert.
5. Discard solution from the lens case immediately and thoroughly clean interior of lens case with a sterile disinfecting solution and air dry after each use.

DESCRIPTION: A sterile, aqueous solution containing proteolytic enzyme (subtilisin) as the active ingredient, and glycerol. Preservative-free.

ACTIONS: This cleaner reduces the adherence of and helps remove protein deposits which can impair the wearing comfort of your lenses

INDICATIONS (Uses): The Boston[®] One Step Liquid Enzymatic Cleaner is indicated for weekly enzymatic cleaning of fluoro silicone acrylate and silicone acrylate rigid gas permeable contact lenses during conditioning (wetting, soaking, and disinfecting), with Boston[®] Conditioning Solution or Boston ADVANCE[®] Comfort Formula Conditioning Solution.

CONTRAINDICATIONS (Reasons not to use): Do not use this solution if you are allergic to any ingredient in this product. Not For Use With Soft (Hydrophilic) Contact Lenses.

WARNINGS:

Do not use directly in the eye. This solution is potentially irritating if misused. Redness or symptoms of irritation may occur if lenses are inadvertently inserted into the eyes without thoroughly rinsing off the enzymatic cleaner. If instilled into the eyes, immediately remove lenses and rinse your eyes with cool fresh water until all of the enzymatic cleaner is thoroughly removed from the eyes. If problems persist, contact your eye care practitioner.

To avoid contamination, never touch dropper tip of the

container with your hands or to any surface.

Contact lenses and solutions, if not handled properly, can become contaminated by common bacteria known to cause ocular infections. Serious damage to the eye and loss of vision may result from using contaminated solutions or from wearing contaminated or damaged contact lenses.

PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYES.

It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use of your lenses and lens care products, including the lens case.

EYE PROBLEMS, INCLUDING CORNEAL ULCER, CAN DEVELOP RAPIDLY AND CAN LEAD TO LOSS OF VISION. THEREFORE, IF YOU EXPERIENCE EYE DISCOMFORT, EXCESSIVE TEARING, VISION CHANGES, OR REDNESS OF THE EYES, IMMEDIATELY REMOVE YOUR LENSES AND PROMPTLY

CONTACT YOUR EYE CARE PRACTITIONER.

Lens Case Care:

- **Clean contact lens cases with digital rubbing with fresh, sterile disinfecting solution. Never use water. Cleaning should be followed by rinsing with fresh, sterile disinfecting solutions (never use water) and wiping the lens cases with a fresh, clean tissue is recommended.** Air-drying or recapping the lens case lids after use without any additional cleaning methods should be avoided. If air drying, be sure that no residual solution remains in the case before allowing it to air dry.
- Replace the lens case at least once every 3 months or as directed by your eye care practitioner.
- Contact lens cases can be a source of microbial growth.

WARNING:

- Do not store your lenses or rinse your lens case with water or any non-sterile solution. Only use fresh solution so you do not contaminate your

lenses or lens case. Use of non-sterile solution can lead to severe infection, vision loss or blindness.

Water Activity:

- Do not expose your contact lenses to water while you are wearing them.
- Ask your eye care practitioner for recommendations about wearing your lenses during any activity involving water.

Important Safety Information:

- Always follow the product directions for use. Failure to follow product directions may lead to vision loss.
- Visit your eye care practitioner regularly.
- Always wash and dry hands before handling lenses.
- Do not use tap water, bottled water or saliva with lenses or lens case.
- Only use fresh solution to clean and disinfect contact lenses.
- Always discard any remaining solution in your lens case after each disinfection cycle.
- Saline solution or rewetting drops will not disinfect your lenses.
- Always replace your solutions, lenses and lens case as directed.
- To avoid contamination, do not touch tip of container to any surface.
- Not for use with heat (thermal) disinfection.

Precautions: To ensure proper hygiene and handling of contact lenses and contact lens solutions, always follow label directions carefully.

- Always discard solution from lens case after each use.
- Discard container immediately after use.
- Store solution at room temperature. Do not freeze.
- Use before the expiration date marked on the carton.
- Keep out of reach of children.
- Never soak lenses in Boston® One Step Liquid Enzymatic Cleaner Only.
- Never place more than 2 drops of Boston® One Step Liquid Enzymatic Cleaner into each lens well.

Adverse Effects: The following problems may occur:

- Eyes stinging, burning, itching (irritation)
- Excessive watering (tearing) of the eyes
- Unusual eye secretions
- Redness of the eyes

- Reduced sharpness of vision (visual acuity)
- Blurred vision
- Sensitivity to light (photophobia)
- Dry eyes

If you notice any of the above conditions, immediately remove and examine lenses. If the problem stops and lenses appear to be undamaged, thoroughly clean, rinse and disinfect the lenses and reinsert.

If the problem continues IMMEDIATELY remove the lens or if a lens appears to be damaged, DO NOT REINSERT. Contact your eye care practitioner immediately.

If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and treatment to avoid serious eye damage.

How supplied: Boston® One Step Liquid Enzymatic Cleaner is supplied in 5 or 15 sterile 0.01 FL OZ (0.3 mL) single-use dispensers.

Use only if tab and single-use container are intact.

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.

© 2021 Bausch & Lomb Incorporated or its affiliates

**Manufactured by: Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA**

Questions or Comments:

1-800-333-4730

www.bausch.com

Made in France