

IMPORTANT: Please read carefully and keep this information for future use.

BAUSCH + LOMB

Boston[®] Rewetting Drops

For Fluoro Silicone Acrylate and Silicone Acrylate Rigid Gas Permeable Contact Lenses

DIRECTIONS FOR USE:

1. Use while contact lenses are on the eyes.
2. Instill 1 to 3 drops into each eye as needed or as directed by your eye care practitioner.
3. Discard solution within 90 days after opening. Record date opened in space provided on bottle label.

DESCRIPTION: Boston[®] Rewetting Drops is a sterile, aqueous, buffered, slightly hypertonic solution containing hydroxyethyl cellulose, polyvinyl alcohol, and a cationic cellulose derivative polymer; preserved with chlorhexidine gluconate (0.006%) and edetate disodium (0.05%).

ACTIONS: Boston Rewetting Drops enhances lens wettability to reduce friction against the cornea. It also helps remove debris that may cause irritation, discomfort and blurring during lens wear.

INDICATIONS (Uses): Boston Rewetting Drops is indicated for use directly in the eye to lubricate and rewet fluoro silicone acrylate and silicone acrylate rigid gas permeable contact lenses.

CONTRAINDICATIONS (Reasons Not to Use): Do not use this solution if you are allergic to any ingredient in this product. **Not For Use With Soft (Hydrophilic) Contact Lenses.**

Contact lenses and solutions, if not handled properly, can become contaminated by common bacteria known to cause ocular infections. Serious damage to the eye and loss of vision may result from using contaminated solutions or from wearing contaminated or damaged contact lenses.

PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYES. It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown that the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement. Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes or redness of the eyes, immediately remove your lenses and promptly contact your eye care practitioner. It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

Important Safety Information:
WARNINGS: To avoid contamination of your solution or your lenses:

- Never touch dropper tip of the container with your hands or to any surface.
- Keep bottle tightly closed when not in use.

PRECAUTIONS: To ensure proper hygiene and handling of contact lenses and contact lens solutions, always follow label directions carefully.

- Always wash your hands before handling your lenses.
- Store solution at room temperature. Avoid freezing.
- Use before expiration date on the container.
- Discard any remaining solution within 90 days after opening.
- Keep out of reach of children.

ADVERSE EFFECTS:

The following problems may occur:

- Eyes stinging, burning, itching (irritation)
- Excessive watering (tearing) of the eyes
- Unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (visual acuity)
- Blurred vision
- Sensitivity to light (photophobia)
- Dry eyes

If you notice any of the above conditions:

IMMEDIATELY remove and examine the lenses.

- If the problem stops and the lenses appear to be undamaged, thoroughly clean, rinse and disinfect the lenses and reinsert.
- If the problem continues or if a lens appears to be damaged, DO NOT REINSERT. Contact your eye care practitioner immediately.

If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and treatment to avoid serious eye damage.

When you visit your eye care practitioner for the assessment of a potential eye infection take your lenses, solutions and lens case to your appointment for possible culturing to better establish the identity of any microorganisms.

GOOD LENS CARE HYGIENE:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens first, the right or the left, to avoid mix-ups.

HOW SUPPLIED: Boston[®] Rewetting Drops is supplied sterile in 0.33 FL OZ (10 ml) bottles. All retail Boston[®] Solutions are packed in unit cartons. Bottles are sealed with an imprinted neckband and marked with a lot number and expiration date.

Questions or Comments:

1-800-553-5340

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.

© 2021 Bausch & Lomb Incorporated or its affiliates

Manufactured by:

Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

www.bausch.com