

**Please read carefully and keep this package insert for future reference.
It is recommended that all eye care professionals review this information with their patients.**



Conditions, Cleans, Removes Protein, Disinfects, Rinses, and Stores.

Directions:

To condition, clean, remove protein and disinfect lenses, complete these simple steps:

- STEP 1:** Place at least 3 drops of **renu[®] Advanced Formula multi-purpose solution** on each side of lens surface and gently rub for 20 seconds.
- STEP 2:** Thoroughly rinse each side of the lens for 5 seconds with **renu Advanced Formula multi-purpose solution**.
- STEP 3:** Place cleaned contact lenses in the lens case and fill with fresh **renu Advanced Formula multi-purpose solution**. Soak at least 4 hours. Remember to always use fresh solution – discard solution from lens case after each use.

Your lenses are now ready to wear. If any debris remains on contact lenses, rinse with **renu Advanced Formula multi-purpose solution** prior to insertion.

Always follow your eye care professional's instructions. Based upon your individual tear chemistry and lens-wearing schedule, your eye care professional may recommend additional products or procedures.

If not wearing contact lenses immediately, store in closed lens case. Do not store your lenses in simple saline in place of **renu Advanced Formula multi-purpose solution**. Saline solution will not disinfect.

Lenses may be stored in the unopened case until ready to wear, up to a maximum of 30 days. If you store your lenses for longer periods of time, they must be cleaned and disinfected with fresh **renu Advanced Formula multi-purpose solution** prior to lens insertion.

Contents:

A sterile, isotonic solution that contains poloxamine, poloxamer 181, diglycine, sodium citrate, boric acid, sodium borate, edetate disodium and sodium chloride and preserved with a triple disinfectant system (polyaminopropyl biguanide 0.00005%, polyquaternium 0.00015% and alexidine 0.0002%).

Actions:

When used daily, **renu Advanced Formula multi-purpose solution** dissolves protein and cleans, loosens and removes accumulations of film, debris and deposits from soft contact lenses. **renu Advanced Formula multi-purpose solution** helps prevent the formation of irritating deposits on the lens surface. It kills harmful microorganisms on the lens. **renu Advanced Formula multi-purpose solution** can also be used to rinse lenses.

renu Advanced Formula multi-purpose solution Benefits:

renu Advanced Formula multi-purpose solution has a unique triple disinfectant system. **renu Advanced Formula multi-purpose solution** rehydrates lenses.

Indications (Uses):

renu Advanced Formula multi-purpose solution is indicated for use in the daily conditioning, cleaning, removal of protein deposits, rinsing, chemical (not heat) disinfection, and storage of soft (hydrophilic) contact lenses, including silicone hydrogel contact lenses, as recommended by your eye care practitioner.

Contraindications (Reasons not to use):

There are no known contraindications for use of **renu Advanced Formula multi-purpose solution**. However, if you are allergic to any ingredient in this product, do not use.

WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN CORNEAL INFECTION AND/OR ULCER AND LEAD TO LOSS OF VISION. It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use of lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care professional. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement. Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care professional. It is recommended that contact lens wearers see their eye care professional twice each year, or if directed, more frequently.

Soaking and Storing Your Lenses:

- You should fill your lens case with fresh solution every time you store your lenses, and never "top-off" or re-use solution. You should discard your solution immediately after your lenses have been removed from the lens case.
- You should not expose or store your lenses in or rinse your lens case with any water, such as tap, bottled or distilled, or with any non-sterile solution.

WARNING:

- Failure to discard solution from lens case after each use, "topping off" solution or use of water to care for your lenses may lead to contamination resulting in eye injury and potential loss of vision.

Rub and Rinse Time:

- Rub your lenses for 20 seconds using at least 3 drops of solution on each side of lens surface.
- Rinse each side of the lens for 5 seconds with solution.

WARNING:

- You should follow the complete recommended lens rubbing and rinsing times in the product labeling to adequately disinfect your lenses and reduce the risk of contact lens contamination. Reduced rubbing or rinsing time may not adequately clean your lenses.
- Never use water, saline solution, or rewetting drops to disinfect your lenses. These solutions will not disinfect your lenses. Not using the recommended disinfectant can lead to severe infection, vision loss or blindness.

Lens Case Care:

- Clean, rinse and air-dry your lens case each time you remove your lenses. In order to permit excess solution to drain, you can flip over your lens case while air drying.
- Replace your lens case frequently, depending upon your hygiene habits.

WARNING:

- Do not store your lenses or rinse your lens case with water or any non-sterile solution. Only use fresh multi-purpose solution so you do not contaminate your lenses or lens case. Use of non-sterile solution can lead to severe infection, vision loss or blindness.

Water Activity:

- Do not expose your contact lenses to water while you are wearing them.
- Ask your eye care professional for recommendations about wearing your lenses during any activity involving water.

Discard Date:

- Discard any remaining solution 90 days after opening.

WARNING:

- Using your multi-purpose solution beyond the discard date could result in contamination of the solution and can lead to severe infection, vision loss or blindness.

Important Safety Information:

- Always follow the product directions for use. Failure to follow product directions may lead to vision loss.
- Visit your eye care practitioner regularly.
- Always wash and dry hands before handling lenses.
- Do not use tap water, bottled water or saliva with lenses or lens case.
- Only use fresh solution to clean and disinfect contact lenses.
- Always discard any remaining solution in your lens case after each disinfection cycle.
- Saline or rewetting drops will not disinfect your lenses.
- Always replace your solution, lenses and lens case as directed.
- To avoid contamination, do not touch tip of container to any surface. Replace cap after using.
- Not for use with heat (thermal) disinfection.

Precautions:

- Always discard solution from lens case after each use.
- Keep the bottle tightly closed when not in use.
- Store solution at room temperature. Do not freeze.
- Use before the expiration date marked on the carton and bottle.
- Discard any remaining solution 90 days after opening.
- Keep out of reach of children.

Adverse Reactions (Problems and what to do): The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

Immediately remove your lenses.

- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, do not put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.
- **If the problem continues, immediately remove the lens and consult your eye care practitioner.**

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

When you visit your eye care professional for the assessment of a potential eye infection take your lenses, solutions and lens case to your appointment for possible culturing in order to better establish the identity of any microorganisms.

Good lens care hygiene:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse the lens case with fresh **renu Advanced Formula multi-purpose solution**, and allow to air dry.

How supplied: **renu Advanced Formula multi-purpose solution** is available for purchase in sterile 2 Fl. Oz. (60 mL), 4 Fl. Oz. (118 mL), 12 Fl. Oz. (355 mL) and 16 Fl. Oz. (473 mL) plastic bottles. Bottles and cartons are marked with a lot number and expiration date.

100% Satisfaction Guaranteed. Your satisfaction with **renu Advanced Formula multi-purpose solution** is guaranteed or your money back. If you are not completely satisfied, need to report an adverse event, or need additional product information call 1-800-553-5340.

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.

©Bausch & Lomb Incorporated

U.S. Patents 8,324,171 and 8,664,180

Manufactured by:
Bausch & Lomb, Rochester, NY 14609