



Let your custom soft potential
FLOURISH

Lenses in small and large diameters— from 10.0 mm to 16.0 mm

Dual Elliptical Stabilization™— an orientation and rotational lens stability technology

Custom base curve range from **6.5 mm to 9.7 mm**

Multifocal and Multifocal Toric Customizable Center Near Zone— **1.8 mm to 3.0 mm**

HOW TO CONFIDENTLY FIT REVIVE™ CUSTOM SOFT CONTACT LENSES

STEP 1

Take Base Measurements to Determine Base Curve

- Horizontal Visible Iris Diameter (HVID)

REVIVE™ BASE CURVE SELECTION CHART													
Average K	50	49	48	47	46	45	44	43	42	41	40	39	38
Lens Diameter	13.5 & Smaller HVID ≤ 11.0 mm	7.7		8.0			8.3			8.6			
	14.0 & 14.5 HVID > 11.0 < 12.5 mm	8.0	8.3			8.6			8.9				
	15.0 & Larger HVID ≥ 12.5 mm	8.6			8.9			9.2			9.5		

STEP 2

Determine Distance Power

- For spherical and toric lenses, use standard vertex adjustment to convert the refractive distance Rx
- For toric lenses, confirm axis orientation

For single vision, order the lenses. To fit multifocal, *continue to step 3*

STEP 3

Determine ADD Power

- For distance powers of plano or minus power, order the refractive ADD power
- For distance powers of +0.25D or greater, order the refractive ADD power less 0.25D
- For distance powers greater than +2.00D and refractive ADD powers greater than +2.25D, order the refractive ADD power less 0.50D Rx

STEP 4

Determine ADD Zone Diameter with Pupil Size

- A standard central ADD zone of 2.4 mm will accommodate most patients
- For pupil sizes smaller than 3.0 mm, a near ADD zone of 1.8 mm to 2.3 mm is recommended
- If pupil size is greater than 4.0 mm, a near ADD zone of 2.5 mm to 3.0 mm is recommended
- If pupils are equal, start with same near zone diameter OU

TROUBLESHOOTING

- If distance vision is blurry, reduce the ADD zone diameter.
- If near vision is blurry, over-refract to obtain the least plus power for acceptable near vision and add this over-refraction to the ADD power. Additionally, the ADD zone diameter can be increased to improve the near vision.
- It is not uncommon to have the ADD zone diameter in the dominant eye set 0.1 mm smaller than the non-dominant eye, even with pupils of equal size.

For more information or questions about custom fitting with Revive™, contact Consultation at 800-253-3669.